



## NEONATOLOGY

Complete range for the smallest patients:

- Ventilators
- Nursing care, thermotherapy and phototherapy
- Accessories

## LIFETHERM 2002/2003

### Resuscitation units

A family of resuscitation units with versatile options for optimal adjustment to the corresponding purpose, the needs of the infant, and the operating staff.

- Radiant warmer with integrated examination lamp
- Conductive heat provided by heated gel mattress<sup>1</sup>
- Quiet electric height adjustment, operation from both sides (up to 30 cm)\*
- Gentle electrical tilting of the bed surface at the touch of a button (Trendelenburg and reverse Trendelenburg positioning)
- Hinged side panels available in different heights
- Skin temperature control and display (accuracy  $\pm 0.1^{\circ}\text{C}$ )\*\*
- Two-sided flip-top or single-piece cover available for bed area\*
- Sturdy rails for attaching trays and ventilator
- Standard rails for mounting accessories such as catheter quiver, storage trays or suctioning accessories\*
- Easy cleaning thanks to optimal design



### ISOTHERM radiant warmer

Universal warmer for thermotherapy, optionally equipped with height-adjustable mobile stand or for ceiling or wall mounting.

- Ceramic, scaling-free heating element with integrated halogen examination lamp (glare-free)
- The reduced surface temperature only has a moderate effect on the room temperature
- Controlled radiation intensity
- Even heat distribution
- Does not emit any visible infrared light (skin colour remains neutral)
- Available in mobile stand design or for wall or ceiling mounting



<sup>1</sup> Gel mattress in conjunction with heated bed surface only

\*optional

\*\* Lifetherm 2002 only



## Sales and Service:

Löwenstein Medical  
Arzbacher Straße 80,  
56130 Bad Ems  
Deutschland  
T: +49 2603 9600-0  
F: +49 2603 9600-50  
info@hul.de  
www.hul.de



To the website:



Images may depict optional accessories that can be ordered separately.

Follow us:



@loewensteinmedical

